

Activity description

Describe your work in a step by step list below, or submit a copy of your synthesis pathway. Please remove the examples before you hand in the form.

If you are doing characterization only, skip this section and fill out section Characterization only.

Operation of every instrument requires a license in LIMS. Licenses are issued after participation on the instrument course in question, where instrument risks are covered. Users must be familiar with the risk matrix for all instruments that they use. Information about all instruments can be found on ntnu.norfab.no.



| Step | Risks involved | Risk reducing measures |
|-----------|---|--|
| 1 | <p>Wafer Clean (if necessary) Piranha solution (3xH₂SO₄ +1x H₂O₂)</p> <p>RCA Clean (if necessary) SC-1: 1xNH₄OH+1x H₂O₂+5x H₂O SC-2: 1x HCl +1x H₂O₂+5x H₂O H₂O rinse and dry</p> | <p>Highly oxidising, reacts violently with organics, highly corrosive, exothermic, can reach 100°C</p> <p>Corrosive, heated to 75°C</p> <p>Fumehood; Wear apron, faceshield and acid gloves; no organics on surface, slowly mix peroxide into acid. Use heavy duty beakers and waste bottle with red cap.</p> <p>PPE as for Piranha. Add acid/base to water, not other way round. Both require buddy</p> |
| 2 | <p>Deposition of 1st metal layer stack Sputter clean 20 nm Cr or Ti adhesion layer 200 nm Pt or Au electrode layer 20 nm Cr or Ti adhesion layer by sputter deposition or e-beam evaporation</p> | |
| 3 | <p>Photolithography Dehydration bake Resist spin (S1823 or AZ9260), softbake at 115°C Exposure Development in MF319 or similar base Descum in O₂ plasma</p> | <p>Solvent exposure</p> <p>Do not inhale fumes</p> |
| 4 | <p>1st metal layer stack etching <u>Dry etch</u>: ICP-RIE using Cl₂ & CF₄ or sputter etching</p> | <p>ICP-RIE: Need to consider platen temperature needs to be above 70°C at least to avoid redeposition of Pt, Au not currently possible – potential issues with resist at higher temperature, may need to use hardmask or sputter etch at lower temperatures instead, or consider wet etch (tried and tested)</p> <p>Keep platen temperature below 70°C (?) if using resist in chamber. Do not exceed 500 V DC bias to avoid burning out Pt thermocouple in chamber, as advised by Oxford instruments</p> |
| <i>or</i> | | |



| | | | |
|----|--|--|--|
| | <p><u>Wet etch:</u> Ti: 1xNH₄OH+2x H₂O₂ (<100 ml) Cr: CR14 (<100 ml) (22% ceric ammonium nitrate, 9% glacial acetic acid, 69% DI water, by weight) Au: 8g KI, 2g I₂, 80 ml DI water or 80 ml of transparent etch (20 g I₂, 80 g NH₄I, 600 ml ethanol, 400 ml dl H₂O) Pt: pre-sputter in Ar plasma 2 min, then Aqua regia at 50°C (25 ml HNO₃+50 ml HCl) Photoresist strip in acetone 10 min O₂ plasma strip</p> | <p>Corrosive Corrosive</p> <p>Highly Corrosive, reacts violently with organics</p> | <p>Apron, Faceshield, acid gloves, extracted wet bench for al the metal wet etches.</p> <p>In addition to above PPE, buddy is presumably required for this, use small amount only. No organics on surface.</p> |
| 5 | <p>O₂ plasma clean PECVD deposition of insulator SiO₂, Si₃N₄ or Al₂O₃</p> | Silane gas | Pump chamber/loadlock three times |
| 6 | Photolithography (S1813) | As before | As before |
| 7 | <p>RIE etch of 1st insulator Using CHF₃, O₂, Ar Photoresist strip in acetone 10 min O₂ plasma strip</p> | Click to write text | Click to write text |
| 8 | <p>Deposition of 2nd metal layer stack Sputter clean 20 nm Cr or Ti adhesion layer 200 nm Pt or Au electrode layer 20 nm Cr or Ti adhesion layer by sputter deposition or e-beam evaporation</p> | Click to write text | Click to write text |
| 9 | Photolithography (S1813) | As before | As before |
| 10 | 2nd metal layer stack etching | As before | As before |
| 11 | 2nd PECVD deposition of insulator layer | As before | As before |
| 12 | <p>NIL PDMS stamp preparation Amonil spin SCIL nanoimprint on MA6 mask aligner</p> | Need additional information on this, I see one of the steps requires an HF etch of the glass carrier plate to remove micro-cracks; is this really necessary? | Work in fume hood, wear suitable gloves and PPE |
| 13 | RIE etch of 2nd insulator layer | As before | As before |
| 14 | <p>2nd metal layer stack nanowell etching ICP-RIE using Cl₂ & CF₄</p> | As before | As before |
| 15 | RIE etch of 1st insulator | As before | As before |
| 16 | Photolithography | As before | As before |
| 17 | Electroplating Ni/Au | Done externally | Done externally |
| 18 | Photolithography | As before | As before |
| 19 | Scribing/Dicing | Small particular matter | Wear gloves and safety specs |

| | | | |
|----|---|---|--|
| 20 | Flip Chip Mounting & Underfill Using solder beads <i>or</i> indium / gallium mix <i>or</i> wirebonder Underfill with UV curable epoxy or Masterbond EP30Med Acetone, IPA, H ₂ O clean O ₂ descum | Gallium – suspected renal toxicity in large doses, potential for contact dermatitis Indium – considered non-toxic in pure form | Very small amounts used |
| 21 | Silanisation Vapour phase silanisation in dry N ₂ or Ar atmosphere using max. 1-2 ml 3-APTES or 3-GEOPS. Bake at 120°C, wash in acetone, IPA, H ₂ O clean O ₂ descum | Silanes: Corrosive, harmful, irritant; Harmful if swallowed. Causes burns. | Safety Goggles, appropriate gloves Use only in extracted wet bench. |

What risks are involved in case power, process gases, N₂, water or pressurized air supply is cut?

Wet etching: Power supply loss would stop extraction of hazardous gases from fume hood.
 Loss of water supply would be hazardous in case of a spill

Equipment: valves actuated through pressurized air would stop to function, loss of water or liquid nitrogen would affect cooling on devices such as the PECVD/RIE and sputter coater.

Process gasses: Loss of nitrogen would prevent correct purge and pumpdown cycle to remove dangerous process gases from chamber and loadlock on PECVD/RIE

Comments to the steps above, if any

| |
|--|
| |
|--|

Characterization

Operation of every instrument requires a license in LIMS. Licenses are issued after participation on the instrument course in question, where instrument risks are covered. Users must be familiar with the risk matrix for all instruments that they use. Information about all instruments can be found on ntnu.norfab.no. Please remove the examples before you hand in the form.

Information about all instruments listed can be found on ntnu.norfab.no

| Instrument | Material to characterize | Comments |
|------------------------------|---|--|
| S(T)EM | Nanoholes through Au or Pt on Si ₃ N ₄ , SiO ₂ , or Al ₂ O ₃ | Sample size can be cut to fit instrument. Conducting and non-magnetic. |
| Dektak or other Profilometer | Photoresist and etched metals/insulators on Si wafers | |
| Optical Microscope | Photoresist and etched metals/insulators on Si wafers | |
| Interferometer | Photoresist and etched metals/insulators on Si wafers | |
| AFM | 200 nm or larger Nanoholes through metal insulator sandwich structure | |



| Comments to the information above, if any |
|---|
| |

Chemicals

Only "Name of chemical" is needed if you do not know the other details. The chemical archive on ecoonline.no shows the chemicals that are provided by NanoLab, and if buddy/logging is needed. For guest log in: Code 1560, username and password "nano". Please remove the examples before you hand in the form.

If the chemicals will be heated, pressurized or mixed with other chemicals, a risk evaluation for this has to be submitted. NTNU's standard risk evaluation forms (found on LIMS, under Info) or any other format of your choice can be used.

NanoLab's licenses related to chemicals (chemical handling license + area courses) do not give any education in safe handling of chemicals. Supervision or previous experience is a prerequisite for chemical work in the cleanroom.



| Name of chemical | Provided by NanoLab | Buddy needed | Logging needed | Will be heated, pressurized or mixed with other chemicals |
|--|---|---|---|---|
| Ethanol (97%) | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No |
| Conc. Nitric acid | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No |
| Conc. Sulphuric acid | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | <input type="checkbox"/> Yes <input type="checkbox"/> No |
| Conc. Hydrochloric acid | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No |
| Ammonium Hydroxide | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No |
| Hydrogen peroxide 30% | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No |
| Tween-20 | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No |
| Biotin-PEG12-NHS Ester | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No |
| Avidin | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No |
| Streptavidin | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No |
| Guanidine hydrochloride | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No |
| Ethylenediaminetetraacetic acid (EDTA) | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No |
| Acetone | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No |
| IPA | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No |
| 3-aminopropyltriethoxysilane | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | <input type="checkbox"/> Yes <input type="checkbox"/> No | <input type="checkbox"/> Yes <input type="checkbox"/> No | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No |
| (3-Glycidyloxy-propyl)trimethoxysilane | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | <input type="checkbox"/> Yes <input type="checkbox"/> No | <input type="checkbox"/> Yes <input type="checkbox"/> No | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No |
| S1813 | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No |
| S1805 | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No |
| S1823 | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | <input type="checkbox"/> Yes <input type="checkbox"/> No | <input type="checkbox"/> Yes <input type="checkbox"/> No | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No |
| AZ9260 | <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No |



Comments to the information above, if any

Silanes will be used in PDMS area
AZ 9260 is introduced in activity NL0132

